

# Sefam S.Box

BY  
STARCK®



A new approach to sleep apnea treatment

T R E A T M E N T





- ▷ Design
- ▷ Ergonomics
- ▷ Technology
- ▷ e-Health

# The S.Box® concept: towards a personalized medicine

Thanks to its multi-connectivity (Bluetooth®, 3G, Wi-Fi) the Sefam S.Box® integrates physiological parameters and data related to the patient's lifestyle in addition to the treatment efficacy indicators.

The S.Box® concept offers healthcare professionals the opportunity to make patient treatment and follow-up more precise as part of an adapted care pathway.



The SEFAM Connect telemonitoring platform assures the security and integration of all information transmitted by the S.Box® and other compatible connected devices.

The mobile apps SEFAM Access (for the patient) and SEFAM Access Pro (for the home care provider) provide access to this information as well as wireless communication with the S.Box®.

The SEFAM Analyze software provides precise and detailed treatment analysis to the healthcare professional.

# S.Box®: an aggregate of technologies and innovations



## iRamp\*

Intelligent ramp detects the moment when the patient is falling asleep



## Intelligent Start

Starts treatment automatically upon breathing in the mask



## Wave & Go

Activates the display upon waving in front of the screen



## Mask Fit & Go\*

Easy check of the mask fit and leakage indication



## CC+ mode\*

Comfort Control Plus  
3 levels of comfort mode



## Circuit Select\*

Automatically calculates the resistance/pressure drop of the entire patient circuit



## Mask Select\*

Calculates precisely the unintentional leak of any mask



## Auto-CPAP mode\*

Reacts to all respiratory events, including flow limitations and snoring



## Adaptive Thermo Control

Automatically adapts the heating of the heated tube to the setting of the heated humidifier



## Humid Control +\*

Enables humidification without condensation



## BT Oximetry

Full SpO2 data acquisition and synchronization with machine



## SEFAM Analyze

New treatment analysis software

(\*) Sefam Patent



# The **patient** at the center of treatment


## Sefam Access mobile app (for patients)

- ▷ Day-to-day follow up of treatment efficacy
- ▷ Integration of other health data from peripheral devices
- ▷ Communication with healthcare professionals
- ▷ Remote control of **S.Box®**


## Sefam Access Pro mobile app (for home care providers)

- ▷ Very close follow up of patient treatment indicators
- ▷ **In app patient population management tools**
- ▷ **Preparation of patient (home) visits**
- ▷ **Access to SEFAM Connect telemonitoring platform**
- ▷ Remote settings of **S.Box®** treatment parameters
- ▷ Communication with the patient

Download the free Sefam Access and Sefam Access Pro apps on your Smartphone or tablet:



The S.Box® design, very low noise level, ergonomics and well-studied packaging make its use and transportation very easy and convenient.



For more information:



### Customer Service

+91 93395 66444  
sales@nidekindia.com

[www.nidekindia.com](http://www.nidekindia.com)

Head Office:  
238B, A.J.C. Bose Road, Kolkata 700 020  
Tel: (033) 4017 7000



BRO155-01 / 08-2017

